

- Fiber Optics should be cleaned at least weekly in a solution of warm water and non abrasive cleaning agent (a small amount of washing up liquid should suffice).

## Cleaning and Maintenance (cont.)

- Soak the optic in the solution and clean with a cloth before rinsing and allow this to dry naturally. Do not place the termination (metal end of the optic) in the solution.
- This may need to be cleaned more frequently depending on the level of usage.

The Fibre Optic Sideglow should be untangled and straightened at least once a week to ensure it is at its optimum condition.

## Fiber Optic Safety

- A risk assessment must be completed for each person who is to have access to the Fiber Optics before use.
- Do not allow the Fibre Optics to be put in a user's mouth under any circumstances.
- Ensure the person using the Fibre Optics is under constant supervision.
- Do not allow the Fibre Optics to be put around a person's neck.

## Troubleshooting

Q. The unit doesn't work.

- A1. Is it switched on? If so,  
 A2. Are the fiber optic strands connected to the lightsource? If so,  
 A3. Has the unit scrolled through its colors? If so,  
 A4. Have you pressed a button on the 8 button wireless controller? If so,  
 A5. Have you recently charged the 8 button wireless controller? so,  
 A6. Please contact us for further assistance.

## Contact Details

For USA please email [info@experia-usa.com](mailto:info@experia-usa.com)

For UK and rest of the world email [support@experia-innovations.co.uk](mailto:support@experia-innovations.co.uk)

Telephone: (800) 882-4045

**Please note:** email is our support team's preferred communication and this takes precedence over telephone calls.

Manufactured proudly in the USA by Experia USA, Inc.


© Copyright 2010-2013 Experia

All rights reserved. Experia is a registered trademark.

**experia** **USA**  
 multisensory for everyone



# Superactive Fiber Optic Bundle



For optimal results, please read and follow instructions carefully  
 Please keep for future reference

**Now Made in the U.S.A. !**



**WARNING: CHOKING HAZARD**

Small Parts: Not suitable for children under 3 years.



**3+ years**

Congratulations on your new purchase! The Superactive LED Fiber Optic Bundle has been designed specifically for people with special needs and should, if used in accordance with these guidelines, offer years of switching learning and fun!

The Superactive LED Fiber Optic Bundle offers opportunities to develop numerous skills including color recognition, cause and effect, and many more!

It is the brightest LED Fiber Optic Lightsource currently available (at time of print) and produces multiple shades of colors when used in calming mode.

The Superactive LED Fiber Optic Bundle comes with its own 8 button controller which allows wireless control from a distance of up to 100 feet, as well as a safe plastic 80" fiber optic harness with 100 strands.

## Features

- Ergonomic design
- Completely silent
- Two operating modes; calming and Interactive
- Latest LED technology. No changing bulbs EVER!
- Qi Charging technology\*\*
- Internal, rechargeable lithium-polymer battery - no batteries to change!
- Cool running
- Hundreds of shades of colors
- Easy maintenance
- Durable & easy to use
- 1 Year manufacturers guarantee (collect and return)

## Pack contents

- 1 x Superactive LED fiber optic light-source
- 1 x 8 button wireless controller
- 1 x Fiber optic sideglow – 2m x 100 strands
- 1 x Low voltage transformer
- 1 x Qi Charging Base\*\*
- 1 x USB Wall Charger

\*\* Qi (pronounced "Chee"), is a wireless charging system

## Getting Started

1. Connect fiber optic 'termination' (black metal part where all the fiber optic strands come together) to the front of the Superactive LED Fiber Optic Lightsource and tighten the bolt to hold the termination in position.
2. Place the Superactive LED Fiber Optic Lightsource on a flat solid floor or mount securely on a wall (**CAUTION:** Fasteners not included. Use appropriate fasteners for your type of wall. Consult a competent contractor if you are unsure).
3. Connect the round XLR connector (silver connector with black cable, right) from the Low Voltage transformer to the socket on the back of the unit.



4. Connect the standard 110V AC plug to a nearby outlet. The Fiber Optics will now illuminate and will scroll through the many colors.
5. We HIGHLY recommend that you fully charge the 8-button controller prior to the initial use (see instructions, below). While each controller is tested, they may be packed for some time in advance of shipping, and will need to be fully charged.
6. Press any of the colored buttons to change the color of the Fiber Optics.

## Modes of Operation

Superactive mode is activated by using the 8 Button Wireless Controller to operate the Fiber Optics.

Calming mode is activated when the 8 Button Wireless Controller has been left dormant for 5 minutes. In this mode the Fiber Optics color sequences automatically.

**Please note: Only use the transformer supplied. Use of a different transformer may result in damage to the product and void the warranty.**

## Charging your Qi-enabled Controller (pronounced "Chee")

1. Your charger is equipped with a built-in wireless charging system and charging base.
2. Simply plug the supplied USB Charger into a wall outlet, and then plug the Qi charging pad into the wall charger.
3. Inside the controller is an inductive coil. This pulls electricity from the charging pad, which charges the battery. This coil is located directly under the RED button of the controller. Simply set the controller down over the charging pad, with the charging pad directly underneath the RED button.
4. There is an LED light right next to the red button, which will illuminate GREEN to indicate that the controller is charging successfully.
5. Simply leave in place to charge. When fully charged, the green LED light will go out.
6. Just remove from the charger to use, and place back on the charger when not in use.

**CAUTION:** The USB wall charger supplied is rated at 2.0 Amps. This is critical for the safe and efficient charging of your wireless controller. If you use a different charger, it **MUST** be at least 2 Amps.

## Cleaning and Maintenance

- The unit should be cleaned with a damp soft cloth. A non abrasive cleaning agent may be used for stubborn stains. Ensure no liquid enters the unit.
- This product has been designed for use by people with special needs, however this does NOT mean it is unbreakable. You must ensure the unit is not hit with anything and that no one pulls or pushes on the unit.